**K-BLEACH 16% HOME**

2 jeringas de peróxido de carbamida 16%  
1 jeringa desensibilizante de nitrato de potasio 5%

**DentalSur2003 S.L.**

Avda. Ortega y Gasset 152  
29006 Málaga (Spain), +34 952 040 502  
hola@ipgdental.com [www.ipgdental.com](http://www.ipgdental.com)

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**984 491 808**

[hola@kiyomidental.com](mailto:hola@kiyomidental.com)

**Rx Only**

Leyenda de símbolos / Meaning of symbols	
	Fabrica   Factory:
	Número de lote   Batch Number
	Fecha de validez   Shelf life
	Límites de temperatura Temperature limits
	No lo utilice si el envase está dañado Do not use if packaging is damaged
	Consérvelo seco   Keep it dry
	Manténgalo lejos de la luz solar Keep protected from sunlight
	Fabricante   Manufacturer
	Consulte las instrucciones de uso Refer to instructions for use
	Packaging reciclable   Recyclable packaging
	Peligro: sensibilización cutánea, irritación de piel y ojos   Danger skin sensitization, skin and eye irritation
	Corrosivo   Corrosive

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### Solo uso profesional

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**K-BLEACH 16%** HOME

2 syringes with 16% carbamide peroxide  
1 desensitizing syringe with 5%  
potassium nitrate

**Dentalsur2003 S.L.**

Avda. Ortega y Gasset 152  
29006 Málaga (Spain), +34 952 040 502  
hola@ipgdental.com [www.ipgdental.com](http://www.ipgdental.com)

**Customer Service**

+34 984 491 808

[hola@kiyomidental.com](mailto:hola@kiyomidental.com)

**Rx Only**

## ES Instructions manual

### Professional use only

Please carefully read all the information in this instruction manual before beginning use.

### Description

This kit contains 2 whitening gel syringes (16% carbamide peroxide) of 3 ml each and one 3 ml syringe of desensitizing gel with 5% potassium nitrate, designed to alleviate dental sensitivity. This product is presented as an effective solution to reduce discomfort associated with sensitive teeth.

### Presentation

Whitening: 2 syringes of 3 ml each and 2 applicator tips.  
Desensitizer: 1 syringe of 3 ml and 1 applicator tip.

### Composition

Whitening: Aqua, carbamide peroxide, carbomer, triethanolamine, hydroxypropyl methyl cellulose, xylitol, methylsaliicylate.  
Desensitizer: Aqua, potassium nitrate, xylitol, triethanolamine, carbopol.

### Product Indication

This kit contains everything needed to perform safe dental whitening on all vital teeth, using the at-home whitening technique under the supervision of a dentist.

### Whitening instructions for use

- The dentist, once the splint is made, has to do the first application in the office, explaining to the patient the use and maintenance of the whitening product.
- The patient, after brushing his teeth, will apply the gel every night with the help of a splint.
- The gel is poured onto the buccal side of the splint, the part in contact with the tooth surface.
- After 30-40 minutes with the splint in the mouth, remove and rinse with plenty of water.
- Repeat the process every day until the product is used up.

In case of tooth sensitivity, use the desensitizing syringe, following the Same application process as with the teeth whitening syringe.

### Precautions

- The use of the product must be indicated and always accompanied by a dentist.
- The oral health must be evaluated prior to treatment. Cracks in the enamel, defective restorations, exposed root surfaces and other factors which can endanger the bleaching treatment must be restored before starting the treatment. The product is contraindicated for users who cannot solve the aforementioned problems.

- Not indicated for bleaching teeth with amelogenesis and dentinogenesis imperfecta, severe fluorosis, intense stains by tetracycline and other anomalies of enamel and dentin which put at risk the vitality of the tooth.
  - As a safety precaution, the product should NOT be used during pregnancy or nursing period and under 18 years old.
  - For diabetic users, it is recommended that the disease is properly controlled before and during the use of the bleaching gel in order to avoid irritant effects of the product. Alternatively, consider the use of an in-office bleaching gel which has no contact to gingival tissue.
  - The product is not recommended for people who have recently had a periodontal surgery or gingivitis or in case of known allergy or chemical sensitivity to peroxides, carbopol, etc.
  - Restorations should be adequately sealed and all exposed sensitive dentin should be covered. If a history of sensitivity exists, treat with appropriate restoration or temporarily with dentin sealant.
  - Tell the patient about the level of results which can be obtained with the bleaching. Inform patients that existing crowns, tooth colored fillings, etc., do not whiten and will have to be changed at the end of the treatment (usually 7 days after the bleaching), teeth which are not vital can be bleached less intensely and that a special bleaching may be required.
  - Try to obtain a photograph before and after the bleaching in order to see the results better. It is also a good idea to bleach one arch at a time. This will allow one arch to be a reference for the other one. Tetracycline stains will require more treatment time than other types of disfiguring stains. Teeth with grey colouring, even though it is not caused by using tetracycline, may also not respond well to bleaching.
  - Occasionally a slight purple discoloration may appear in the tray following bleaching adjacent to amalgam fillings.
  - Teach the user how to place the gel in the tray. We recommend one drop of gel in each tooth space of the tray. The excess of gel can be avoided and removed by using a toothbrush or finger. Prolonged contact with the gingiva can result in localized irritation.
  - Provide your patient the leaflet with the instructions and emphasize the importance of reading it.
  - Tell the user about possible sensitivity to temperatura variances, mainly to cold.
  - Exposed root surfaces may experience sensitivity or another adverse factor. Protect the area or cut the tray to the cement union in order to avoid prolonged contact of the gel with the sensitive area. In these cases one cannot expect the uncovered area to bleach with the same intensity. There are cases where the sensitivity is not easy to treat. If none of the treatment forms described have an effect, using the gel can be prescribed for alternating days, every 3 days or for 1 to 2 hours a day prolonging the bleaching. Generally the sensitivity disappears in at most 3 days. If the user has too much discomfort or an allergic reaction to the product, the treatment must be suspended.
  - Guide the user about not ingesting acid drinks or
- not indicated for bleaching teeth with the bleaching. Strongly coloured drinks and food should also be avoided to not affect the bleaching process. Try to follow up the bleaching by asking the patient to return every 5 to 7 days, for example, in order to evaluate the treatment, mainly at the start of the bleaching.
  - Gingival irritation can also occur when using mouthwash or even toothpaste with a strong taste before using the tray with the gel. If you wish to use these products, it is advisable to leave an interval of at least 30 minutes (or sufficient time for the sensitivity to disappear) between using the mouthwash or toothpaste and the gel. Aggressive brushing on the gingiva can also make it sensitive and result in irritation caused by the gel. Gingival irritation can also be caused by a custom tray which was not adequately trimmed. Sharp edges must be trimmed off the tray since it could cause gingival tissue lacerations and affect the safe use of the product.
- ### Side effects
- Occasionally sensitivity may require treatment. If the sensitivity becomes uncomfortable, treat the teeth with a soothing desensitising gel. The more concentration of Peroxide, greater the possibility of suffering dental sensitivity.
  - Areas of hypocalcification can exist which may not be visible to the naked eye. Hypocalcified areas will whiten faster, thereby becoming more obvious during bleaching. In some cases these stains disappear owing to the dehydration of the enamel after concluding the bleaching. Treating these areas with fluoride can also help to make them disappear.
  - Some users can feel slight irritation in the gingiva, tongue or lips, generally caused by excessive use of gel in the tray.
- ### Desensitizer instructions for use
- Brush teeth: Clean your teeth with your regular toothpaste.
  - Apply desensitizing gel to the dental tray: Remove the syringe cap and apply a small amount of gel to the provided dental tray.
  - Keep the tray in the mouth for 30-40 minutes: Gently place the dental tray in your mouth, ensuring the gel is in contact with the sensitive teeth.
  - Remove the tray and rinse: After the recommended time, remove the tray and rinse your mouth with warm water to eliminate any gel residue.
- Repeat this procedure daily until the syringe is completely used. If sensitivity persists, consult your dentist.
- ### Preservation and Storage
- Store the product at 5°C and 25°C / 41°F and 77°F. Do not freeze the product. Keep it out of direct sunlight.
- ### Warning
- Do not use the product after expiry date. Maintain it out of reach of children. The discard of product follows the directive 94/62/CE.
- Irritant to eyes and skin for concentrations over 16% of carbamide peroxide.  
This material was only manufactured for dental use and must be handled according to the instructions for use. The manufacturer is not responsible for damage caused by other uses or incorrect handling. Besides this, the user is obliged to prove, before use and under his or her responsibility that this material is compatible with the use desired, mainly when this use is not indicated in the instructions for use. Data descriptions do not constitute any type of guarantee and so are not in way linked.

Leyenda de símbolos / Meaning of symbols	
	Fabrica   Factory:
	Número de lote   Batch Number
	Fecha de validez   Shelf life
	Límites de temperatura Temperature limits
	No lo utilice si el envase está dañado Do not use if packaging is damaged
	Consérvelo seco   Keep it dry
	Manténgalo lejos de la luz solar Keep protected from sunlight
	Fabricante   Manufacturer
	Consulte las instrucciones de uso Refer to instructions for use
	Packaging reciclable   Recyclable packaging
	Peligro: sensibilización cutánea, irritación de piel y ojos   Danger skin sensitization, skin and eye irritation
	Corrosivo   Corrosive